

Mass Spectrometric Analyses for anti-FLAG Immuno-Precipitated Proteins

2015. 8.

한국표준과학연구원 이진환

Wear powder-free nitrile gloves and lab coat, and perform as many manipulations as possible under the hood.

1. Rinse 100 ml Pyrex bottle with mass spec-grad H₂O.
2. Make up (50 ml for four pull-downs) Lysis buffer in 100 mL Pyrex bottle.
3. Take 50 mL buffer (50 mM HEPES, pH 7.5; 70 mM KOAc; 5 mM Mg(OAc)₂).
4. Add 0.1 g n-dodecyl-D-Maltoside (Thermo Scientific; Cat # 89903).
5. Add 1 Protease Tablet (Roche; Cat # 1117700).

■ **Harvest** (should start with one 90–100% confluent 14 cm dish).

1. Wash cells with in 5 ml Hank's dissociation buffer (Gibco; Cat #13150-016).
2. Add 2 ml Hank's dissociation buffer and incubate in 37 °C incubator for 5 min.
3. Detach cells from the dish by vigorous tapping.
4. Add 13 ml cold PBS and transfer cells to a 50 ml conical tube.
5. Centrifuge.
6. Remove supernatant and add 1 ml cold PBS.
7. Gently resuspend cells and transfer to non-stick microcentrifuge tube (VWR; cat # 20170-650).
8. Spin down and aspirate off media using a syringe needle.
9. Store cells in –80 °C freezer should you want to stop.

■ **Lysis and FLAG Immuno-Precipitation**

1. Add 1 mL lysis buffer to the cell pellet and resuspend cells.
2. Incubate on rotator in cold room for 30 min.
3. Spin in cold microcentrifuge for 15 min at 16,600×g.
4. Prepare beads (anti-FLAG M2 affinity gel; sigma A2220).
 - ① 50 µL beads per sample (~5 mg) in each tube.
 - ② Wash 3 times with 1 ml Lysis buffer.
5. Save 50 µl of the supernatant for Western blot analysis and protein quantitation.
 - For Protein quantification
 - ① Take 5 µl and 95 µL water.
 - ② Add 20 µl to Quick Bradford Buffer.

- ③ Measure OD.
- ④ Calculate concentration—take 30 µg for western blotting and freeze.
6. Add the supernatant corresponding to 5 mg of protein plus lysis buffer to 1 ml to tubes containing beads.
 - ① Incubate 1 hr on rotator at 4 °C.
 - ② Wash 5 times with 1 ml Lysis buffer.
 - ③ Wash 2 times with 1 ml 0.1 M Tris-HCl, pH 8.5 (Mass Spec Buffer).
 - ④ Spin and aspirate off any residual buffer.
 - ⑤ Resuspend in 100 µl MS Buffer.
 - ⑥ Save 5 µl for western blotting.
 - ⑦ Spin and aspirate.

■ Digest

1. Make up 10 M urea.
2. Elute from the beads with 64 µl, place on Thermomixer for 15 minutes.
3. Add to Bio-rad micro bio-spin chromatography column (cat # 732-6204) place in microcentrifuge tube.
4. Spin 2 min at 6,000 X g in microcentrifuge.
5. Add 16 µl MS buffer to dilute to 8 M Urea.
6. Add 0.5 µl 500 mM TCEP.
7. Incubate for 20 min at room temperature.
8. Make up 500 mM 2-chloro-acetamide (18.7 mg in 400 µl MS buffer).
9. Add 1.8 µl.
10. Incubate for 15 min at room temperature.
11. Add 2 µl Lys-C endoprotease (100 ng/ µl).
12. Incubate for 4 hr at 37 °C.
13. Add 240 µl MS Buffer to dilute to 2 M urea.
14. Add 3.2 µl 100 mM CaCl₂.
15. Add 2 µl Trypsin (100 ng/ µl) (Promega; VS11A).
16. Incubate for up to 14 hr at 37 °C and desalt using StageTips (Thermo; SP301) for LC injection.